Splinting materials are intended to be used for fabrication of custom-molded rigid splints, orthoses, and adaptive equipment.

## PRECAUTIONS FOR FINISHED ORTHOSES

**Splint adjustments are to be made only by a healthcare professional. The healthcare professional is responsible for providing wearing instructions and precautions to other healthcare practitioners, care providers involved in the patient's care, and the patient.**

**If unusual swelling, skin discoloration or discomfort occurs, use should be discontinued and a healthcare professional consulted.**

## PATTERNING

Patterns can be drawn on all splinting materials with a scratch awl before placing the material in the water

## CUTTING

All splinting materials can be cut at room temperature by using heavy-duty shears, or scoring with a utility knife and bending at the scored mark, then completing the cut with a knife or scissors. If more detailed pattern-making is required, it is recommended that the pattern be transferred to the splinting material to make a rough blank. Then heat the rough blank for one-half the recommended time or until material is soft enough to cut with scissors. This method will reduce rough edges and save time with edge finishing. Light splinting materials can be cut at room temperature using scissors.

## HEATING

**For a chart giving heating information for each material, see "Water Temperature and Time Guide for Rolyan and Rolyan Aquaplast Splinting Materials"**

The best method for softening all splinting materials is heated water. Use a heat gun only for spot-heating and making adjustments.

- For hot water heating, fill frying pan or heat pan with water at least 1" (2.5 cm) deep. Heat water to temperature recommended in the chart (use thermo meter to measure water temperature). For Aquaplast splinting materials, use a heat pan liner in the bottom of the heat pan or frying pan. HydrocollatorG units can be used for large pieces of splinting material. Use a heat pan liner to dip the material into the unit, to ensure that the material does not touch the bottom. Do not leave the splint blank unattended.
- In addition to hot water heating, San-Splint, Orthoplast, OrthoPlast II and Aquaplast materials can be heated in a convection oven. Note that dry heating makes the material self-bonding, so no surface preparation is required for bonding. Material that is heated in the oven also becomes hotter than material heated in water; check temperature carefully before applying it to the patient. Heat material without

overlapping or folding on a baking sheet with a non-stick surface in a preheated oven as follows.

- San-Splint: heat at 175°F (80°C) for 3 to 4 minutes.
- Orthoplast: heat at 160° to 170°F (70-75°C) for 2 minutes.
- Orthoplast II: heat at 140° to 158°F (65-70°C) for 1 to 2 minutes.
- Aquaplast: heat at 195° to 205°F (90-95°C) for up to 10 minutes or until transparent or translucent.

*Note: Overheating splinting materials increases the draping/stretching characteristics; cool slightly before handling it to avoid excess stretching.*

## FORMING

After the pattern is made, the splinting material should be heated in water as indicated in the temperature and time chart. The entire pattern or piece of material should feel soft when removed from the water. When Aquaplast materials become transparent or translucent in a heat pan, the material is soft and ready to form. When the material is soft, remove it from the water and pat it dry with an absorbent towel. (To avoid towel impressions or lint, cover the towel with a pillow case before drying the splinting material.) **Check the temperature of the material before beginning to fabricate the splint on the patient.**

Forming some splints on a patient is easier with the assistance of gravity. When forming a volar forearm splint, for example, place the patient's hand in supination whenever possible, to allow the splinting material to drape into the arches of the hand.

*Technical Tip: To prevent hands from sticking to Original Aquaplast and Original Resilient, use cold cream or mineral oil on hands or dip hands in cold water.*

Use smooth strokes when molding splints. Original Aquaplast®, Aquaplast-T, Watercolors®, ProDrape® -T, Polyform®, Polyflex II®, Orthoplast II® and TailorSplint® require minimum shaping with gentle strokes. Repeated pressure is not necessary to retain form while cooling. The controlled stretch of Original Resilient®, Resilient-T, Ezeform®, Synergy®, San-Splint® and Orthoplast® allows more aggressive handling. These materials require firm pressure to form the splint. All Aquaplast splinting materials have 100% elastic memory to allow repeated molding for easy and economical splint revisions. Re-soften the material and it will return to its original size and shape.

## COOLING

Cold water, cold packs, cold spray, or cold wraps can be used to speed the cooling process.

*Technical Tip: Exercise band that has been stored in a freezer can also be used to speed the cooling process.*

## ADJUSTMENTS

Change in contour or design of the splint can be accomplished easily by spot-heating the area by pouring

hot water over it or using a heat gun. Care should be taken not to overheat the material when using a heat gun—keep the hot air stream of the heat gun moving and do not hold it too close to the material.

## ADHESIVE PROPERTIES

Aquaplast-T, Resilient-T, ProDrape-T, Watercolors, Polyform, Polyflex II, Orthoplast, OrthoPlast II and TailorSplint have a protective coating to prevent accidental self-bonding. To bond these materials, remove the coating by scraping the surface or brushing the surface with an acetone material. Spot-heat both surfaces with a heat gun until they become tacky, make sure both surfaces are dry, and then press firmly together. This technique will provide a permanent bond, ideal for adding outriggers and other attachments to the finished splint.

San-Splint resist bonding, and using a solvent to prepare the surface for permanent bonding is recommended.

Original Aquaplast, Original Resilient, Ezeform and Synergy do not have a protective coating, and they are self-bonding when both pieces are hot and dry. A solvent is not required—but cleaning the splint surface prior to bonding is recommended.

## Edge finishing

Smooth edges by heating the rough edges of the splinting material with a heat gun or dipping the edge into heated water for 5 to 10 seconds. Smooth the heated edges using finger tips or the palm of the hand—a small amount of hand or body lotion on the finger tips makes smoothing the edges easier.

## CARE AND CLEANING

### Splinting materials

The shelf life of low-temperature splinting materials is indefinite under the following conditions:

- Store at temperatures between 40° and 90°F (4° and 32°C).
- Store at less than 65% relative humidity.
- Avoid prolonged exposure to light, especially ultraviolet light.
- Avoid exposure to corrosive fumes and ethylene oxide fumes.

### Fabricated splints

The splint will lose its shape in temperatures over 135°F (57°C). It should be kept away from sources of heat such as ovens, hot water, open flames, and sunny car windows.

The splint can be cleaned with soap and lukewarm water. The straps can be washed with soap and water, but the adhesive portion should not be moistened. The splint and straps should be allowed to dry thoroughly before reapplication of the splint.

**Please Note:** Characteristic markings on most sheet material include gate and pin marks. The markings are part of the normal production process, are within specification, and are not considered defects. The gate is a circular mark usually in the center of the sheet but may appear as a half moon shape at the edge of the sheet. Pin marks are left behind when the sheet material is pushed off of the mold.

Serious events experienced from using this product should be reported immediately to Performance Health and to the local Competent Authority.

Symbols Glossary located at <http://www.performancehealth.com/medical-device-glossary>.

These user instructions are available in large print, please contact Performance Health Customer Services for a copy.

Please retain these instructions for future reference.

## WATER TEMPERATURE AND TIME GUIDE FOR ROLYAN AND ROLYAN AQUAPLAST SPLINTING MATERIALS

| Water Temperature <sup>1</sup>                      | Approximate Heating Time |                         |                       |                     | Approximate Working Time |                       |                       |                     |
|---|--------------------------|-------------------------|-----------------------|---------------------|--------------------------|-----------------------|-----------------------|---------------------|
|   | 1/8" (1.6 mm)            | 3/32" (2.4 mm)          | 1/4" (3.2 mm)         | 3/16" (4.8 mm)      | 1/8" (1.6 mm)            | 3/32" (2.4 mm)        | 1/4" (3.2 mm)         | 3/16" (4.8 mm)      |
| Polyform 150–160°F (65–70°C)                        | 30 sec. <sup>4</sup>     | —                       | 1 min. <sup>5</sup>   | —                   | 1 min. <sup>5</sup>      | —                     | 3–5 min. <sup>5</sup> | —                   |
| Ezeform 160–170°F (70–75°C)                         | —                        | —                       | 1 min. <sup>5</sup>   | —                   | —                        | —                     | 4–5 min. <sup>5</sup> | —                   |
| Orthoplast II 140–158°F (65–70°C)                   | —                        | —                       | 1 min. <sup>5</sup>   | —                   | —                        | —                     | 3–5 min. <sup>5</sup> | —                   |
| Polyflex II 150–160°F (65–70°C)                     | 30 sec. <sup>4</sup>     | 30–45 sec. <sup>4</sup> | 1 min. <sup>5</sup>   | —                   | 1 min. <sup>5</sup>      | 2–3 min. <sup>5</sup> | 3–5 min. <sup>5</sup> | —                   |
| TailorSplint 150–160°F (65–70°C)                    | 30 sec. <sup>4</sup>     | 30–45 sec. <sup>4</sup> | 1 min. <sup>5</sup>   | —                   | 1 min. <sup>5</sup>      | 2–3 min. <sup>5</sup> | 3–5 min. <sup>5</sup> | —                   |
| Synergy 160–170°F (70–75°C)                         | —                        | —                       | 1 min. <sup>5</sup>   | —                   | —                        | —                     | 4–6 min. <sup>5</sup> | —                   |
| San-Splint <sup>2</sup> 160–175°F (70–80°C)         | —                        | —                       | 90 sec. <sup>4</sup>  | —                   | —                        | —                     | 4–5 min. <sup>5</sup> | —                   |
| Orthoplast 160–170°F (65–75°C)                      | —                        | —                       | 90 sec. <sup>4</sup>  | —                   | —                        | —                     | 4–5 min. <sup>5</sup> | —                   |
| Aquaplast-T <sup>3</sup> 160–170°F (70–75°C)        | 35 sec. <sup>4</sup>     | 60–70 sec. <sup>4</sup> | 1–2 min. <sup>5</sup> | 2 min. <sup>5</sup> | 2 min. <sup>5</sup>      | 3 min. <sup>5</sup>   | 4 min. <sup>5</sup>   | 7 min. <sup>5</sup> |
| Original Aquaplast <sup>3</sup> 160–170°F (70–75°C) | 35 sec. <sup>4</sup>     | 60–70 sec. <sup>4</sup> | 1–2 min. <sup>5</sup> | 2 min. <sup>5</sup> | 2 min. <sup>5</sup>      | 3 min. <sup>5</sup>   | 4 min. <sup>5</sup>   | 7 min. <sup>5</sup> |
| Resilient-T <sup>3</sup> 160–170°F (70–75°C)        | 35 sec. <sup>4</sup>     | —                       | 1–2 min. <sup>5</sup> | —                   | 2 min. <sup>5</sup>      | —                     | 4 min. <sup>5</sup>   | —                   |
| Original Resilient <sup>3</sup> 160–170°F (70–75°C) | 35 sec. <sup>4</sup>     | —                       | 1–2 min. <sup>5</sup> | —                   | 2 min. <sup>5</sup>      | —                     | 4 min. <sup>5</sup>   | —                   |
| ProDrape-T <sup>3</sup> 160–170°F (70–75°C)         | —                        | 60–70 sec. <sup>4</sup> | 1–2 min. <sup>5</sup> | —                   | —                        | 3 min. <sup>5</sup>   | 4 min. <sup>5</sup>   | —                   |
| Watercolors <sup>3</sup> 160–170°F (70–75°C)        | 35 sec. <sup>4</sup>     | 60–70 sec. <sup>4</sup> | 1–2 min. <sup>5</sup> | —                   | 2 min. <sup>5</sup>      | 3 min. <sup>5</sup>   | 4 min. <sup>5</sup>   | —                   |

<sup>1</sup> Below 150°F (65°C), Rolyan splinting materials will not soften sufficiently to become workable.

<sup>2</sup> San-Splint Orthoplast & Orthoplast II can be heated in a convection oven (see "Heating").

<sup>3</sup> Aquaplast materials can be heated in a convection oven (see "Heating"). Heating and working times of Aquaplast materials determined using a water temperature of 160°F (70°C).

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 Performance Health Supply, LLC  
W68 N158 Evergreen Blvd  
Cedarburg, WI 53012, USA  
+1800-323-5547  
www.performancehealth.com

**Distributed in the United Kingdom By:**  
Performance Health International LTD  
Nunn Brook Road  
Huthwaite, Sutton-in-Ashfield  
Nottinghamshire, NG17 2HU, UK  
+44 (0) 1623 448 706  
www.performancehealth.co.uk

**Distribué en France par:**  
Performance Health France  
13 rue André Pingat CS 10045  
51 724 Reims Cedex France  
00(33) 03 10 00 79 30  
www.performancehealth.fr

 **EC REP** CEpartner4U  
Esdoornlaan 13  
3951 DB Maarn  
The Netherlands  
www.cepartner4u.com

 Performance Health France  
13 rue André Pingat CS 10045  
51 724 Reims Cedex France

**Distributed in Canada By:**  
Performance Health Canada Inc  
6715 Millcreek Drive, Unit 2  
Mississauga, ON L5N 5V2 Canada  
+1800-665-9200  
www.performancehealth.ca

**UK Responsible Person:**  
UKCApartner4U Ltd  
7 Campion Way  
Bingham  
Nottingham  
NG13 8TR  
UK  
www.UKCApartner4U.com

 **CH REP** MedEnvoy Switzerland  
Gotthardstrasse 28  
6302 Zug  
Switzerland  
www.medenvoyglobal.com