

# ASPEN SUMMIT™ 456 TLSO ADJUSTABLE

THORACIC LUMBAR SACRAL ORTHOSIS  
LUMBAL-SAKRAL-THORAKALE ORTHESE  
ÓRTESES TORÁCICA LUMBOSACRA  
ORTESI TORACICA SACRO-LOMBARE  
ORTHÈSE THORACIQUE LOMBO-SACRÉE  
LUMBALE SAKRAL THORAKALE ORTHÈSE  
BEDERNÍ KRÍŽOVÁ HRUDNÍ ORTÉZA  
ΛΕΝΔΗΝΙΡVΛΕΣΩΛΕΝ THORAXORTOSE  
ΟΣΦΥΪΚΗ ΘΩΡΑΚΙΚΗ ΟΡΘΩΣΗ ΙΕΡΟΥ ΟΣΤΟΥ  
LUMBAL BRYSTORTOSE  
ORTEZA PIERSIOWO-LĘDŹYWIOWO-KRZYŻOWA  
LANNERANGAN SAKRALAINEN RINTAKEHOORTOOSI  
LUMBAL SAKRAL THORAKAL ORTOS  
ОРТЕЗ ДЛЯ ГРУДНОГО И ПОСНИЧНО-КРЕСТЦОВОГО ОТДЕЛОВ  
ПОЗВОНОЧНИКА  
ÓRTESE TORÁCICA SACRO LOMBAR  
腰骶胸部矫形器

1A 1B  
1C  
2  
3A 3B 3C  
4A 4B 4C 4D  
5

4A 4B 4C 4D  
5

4A 4B 4C 4D  
5

4A 4B 4C 4D  
5

4A 4B 4C 4D  
5

4A 4B 4C 4D  
5

**ENGLISH**  
**1. BELT SIZING**  
A. Fit male tails of belt through proper sizing adjustment holes in sides of back panel.  
B. For smaller sizes, slide side panel through sizing adjustment slot and then fit tabs through proper adjustment holes from the under side.  
C. Once proper size is achieved, pull taught to lock tabs in place.

**2. BACK PANEL**  
Custom form the malleable aluminum support in the back panel to properly fit patient anatomy.

**3. PLACEMENT**  
A. Apply the fully stretched brace making sure the back panel is centered and low on the back.  
B. Overlap the ends of the brace so the hook material grips the opposite side loop as tightly as possible.  
C. Ends may be angled to accommodate for different anatomies.

**4. TIGHTEN AND ADJUST**  
A. Adjust the vertical height of the posterior adjustment strap so that the yoke rests at C7, or 2" below the shirt collar.  
B. Slip each thumb into the pull tabs. Extend the tabs out and away from the back. When desired compression is achieved, place the hook side of the pull tabs down onto the loop material to secure.  
C. Pull and tighten the shoulder straps so they fit wide on the shoulders.  
D. Fasten and tighten the chest buckles to secure the shoulder straps.

**5. MODIFICATION OPTIONS**  
This product may be substantially modified by an individual with expertise. For more information, please visit [www.aspenmp.com/education/product-modifications](http://www.aspenmp.com/education/product-modifications).

**INTENDED USE**  
To provide motion restriction and stability of the Thoracic-lumbar-sacral spine.

**INDICATIONS**  
Postoperative support, Degenerative disc disease, Bulging or herniated disc, Fracture management, Spondylolysis, Spondylolisthesis, Spinal stenosis, Kyphosis, Spondylolysis (Osteoarthritis/Degenerative Joint Disease).

**CONTRAINDICATIONS**  
Hypersensitivity or allergies to any of the materials from which the brace is made. Swelling of the lymphatic tissue caused by circulatory disorders. Patients who are not allowed to wear a brace according to medical instructions.

**PRECAUTIONS**  
READ INSTRUCTIONS BEFORE USE. Proper training in the use of this device should take place before it is applied. These directions are guidelines only and are not offered as medical recommendations. If you suffer from a serious medical condition, we strongly suggest that you consult with a licensed health care professional before using this product. Proper fitting is required for this product to be effective. Under some circumstances, this product may be prescribed by a physician. Please see the limited warranty for further information.

**CARE INSTRUCTIONS**  
Hand wash only at 30° C. Wash with mild detergent; Air dry. Do not use bleach or other chemicals.

**MATERIALS**  
Nylon, Polyester, HDPE, PVC, Polyurethane, UHMWPE.

**NOT MADE WITH NATURAL RUBBER LATEX.**

**DISPOSAL**  
Product may be safely disposed in accordance with local laws.

**REPORTING**  
Provide report any serious incident involving the use of this device to both the Competent Authority in your state and to the manufacturer (listed on this document).

**LIMITED WARRANTY**  
Aspen Medical Products, LLC, Irvine, CA 92618, warrants to the user who originally purchases this product that it is free from defects in material and workmanship. The sole obligation of Aspen Medical Products, LLC in the event of breach of warranty shall be to repair or replace the defective product or parts).

Aspen Medical Products, LLC shall have no obligation under this limited warranty in the event:  
(a) The product was not purchased from Aspen Medical Products, LLC or through its authorized channels of distribution;  
(b) The product is altered;  
(c) Any parts not supplied by Aspen Medical Products, LLC are inserted into the product;  
(d) The product is not used in accordance with the Aspen Medical Products, LLC Instructions for Use.

THE FOREGOING IS THE SOLE AND EXCLUSIVE REMEDY FOR FAILURE IN SERVICE OR, OR DEFECT IN THE PRODUCT. ASPEN MEDICAL PRODUCTS SHALL NOT BE LIABLE UNDER THIS OR ANY IMPLIED WARRANTY FOR ANY DIRECT, SPECIAL, INCIDENTAL, OR CONSEQUENTIAL DAMAGES. THIS WARRANTY IS IN LIEU OF ALL OTHER WARRANTIES, EXPRESSED OR IMPLIED, INCLUDING THE WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE OR USE. AND ALL OBLIGATIONS OR LIABILITIES ON THE PART OF ASPEN MEDICAL PRODUCTS FOR DAMAGES ARISING OUT OF OR IN CONNECTION WITH THE USE OF THE PRODUCTS, WHICH ARE HEREBY DISCLAIMED AND EXCLUDED BY ASPEN MEDICAL PRODUCTS, LLC.

This warranty gives you specific legal rights and you may have other rights which vary from location to location.

A listing of authorized distributors in your area is available upon request.

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Společnost Aspen Medical Products, LLC nebude mít podle této omezené záruky žádnou povinnost v případě, že:
(c) Do produktu jsou vložené jakékoli součásti, které nedodává společnost Aspen Medical Products, LLC; nebo
(d) Produkt se nepoužívá v souladu s pokyny k použití společnosti Aspen Medical Products, LLC.